



DatPass v4.0 Installation CD Read Me

March 18, 2012

fourtec - Fourier Technologies

Please read the following important guidelines relating to installation of the DatPass audit trail database.

- [First-time DatPass Users](#)
- [Users with Existing DatPass Installation](#)
- [Uninstalling Existing DatPass v2.0](#)
- [Move Database Password](#)

First-time DatPass Users

When installing DatPass on your computer for the first time (i.e. no previous DatPass installation exists), the DatPass database (file name: Current.mdb) will be created in the following default path:

\Users\Public\fourtec\DatPass

Users with DatPass Admin privileges can move the database to another location via the DatPass Move Database feature.

Windows Vista Users:

When launching DatPass for the first time, if the database is not found by the software it is because it has not been created in the default path during DatPass installation. The user will then be prompted to locate the missing database, where a copy can always be found in the DatPass directory on \Program Files\fourtec\DatPass.

It is recommended to move the database to either the default path of \Users\Public\fourtec\DatPass or to another location of your choice.

Users with Existing DatPass Installation

Existing users of DatPass v2.0 and 3.0 are able to upgrade to DatPass v4.0 and continue to use their existing audit trail database.

During installation of DatPass v4.0, the software will detect previous DatPass installation on the computer. A new database will therefore not be created. However, DatPass v4.0, when first launched, will attempt to open the existing database from the default path of \Users\Public\fourtec\DatPass. As the existing database is not located here, the user will then



be prompted to locate the missing database. The user simply needs to browse to the path of the existing database, and DatPass v4.0 will then launch with the existing database and audit trail intact.

Uninstalling Existing DatPass v2.0

It is not recommended to maintain more than one installation of DatPass on the computer. Therefore, previous versions should be uninstalled only after DatPass v4.0 is successfully installed.

Note that uninstalling DatPass v2.0 will also delete the Current.mdb file from the v2.0 default database path of My Documents/Fourier Systems/DatPass.

For any installation where the database is located in this default path, it is critical that prior to commencing the uninstall process, the user backs up the Current.mdb or simply moves the database to another location, to avoid having it deleted following uninstall.

Move Database Password

In DatPass v4.0 the password has been changed to **fourtecfr**

Please consult the DatPass User Guide for further software installation and activation instructions.

For technical support please contact support@fourtec.com.