



## **fourtec Declaration of Accuracy**

### *For PicoLite-16K Temperature Data Logger*

The PicoLite data logger is produced using Surface Mounted Device (SMD) assembling technique. Each product has been through complete quality assurance and calibration testing measures. All Fourtec production procedures are fully ISO 9001:2008 certificated.

The PicoLite-16K Temperature USB data logger is rated to the following temperature specifications:

Range: -40 to 40 °C

Accuracy:  $\pm 0.5$  °C across the range

For further information please contact [support@fourtec.com](mailto:support@fourtec.com)

fourtec – Fourier Technologies.  
Quality Assurance Department



**Manufactured under ISO 9001:2008**

[www.fourtec.com](http://www.fourtec.com)