



## **DatPass Admin Software Release Notes**

**FDA 21 CFR Part 11 Compliant**

### ***System Requirements***

#### ***Software***

- Windows XP SP3, Windows Vista, Windows 7 (32-bit and 64-bit)
- Minimum screen resolution: 1024 x 768 (800 x 600 not supported)

#### ***Hardware***

- Pentium 4, 2.0 GHz or better
- 512 MB RAM
- 250 MB available disk space for the DatPass application

**Release Date:** 18 March, 2012

**DatPass Software Version:** 4.0.0.0

**Software Support:** DataSuite, DataNet CFR, MicroLab CFR

**Hardware Support:** DataNet, DaqLink, MicroLogPRO II (Rev A and Rev B), MicroLite I & II

### ***Release Information***

- Added support for DataSuite CFR software
- Software access using USB security dongle no longer required.
- Software activation mechanism added, involving free 30-day trial and activation with a unique serial number license code (SN format: 9901010XXXXXX)
- Added ability to copy and paste data direct from the Audit Trail to Windows® clipboard.
- Added ability to filter Audit Trail column headers
- DatPass User Guide available from Help menu

### ***Contacting fourtec Technical Support***

**e-mail:** [support@fourtec.com](mailto:support@fourtec.com)

**Web:** [www.fourtec.com](http://www.fourtec.com)

**Voice:** International +972-3-901-4849

Within USA only (toll-free) 1-888-778-4927