



## MicroLab and MicroLab Lite Validation Certificate

This document certifies that the Fourier Systems software MicroLab, supporting the MicroLog, MicroLogPro and MicroLite data loggers and MicroLab Lite, supporting MicroLite, have been validated to display correct data. Both MicroLab and MicroLab Lite display data downloaded from the aforementioned data loggers via USB or Serial cable connection to the PC. The internal or external sensor data is displayed in graph, tabular and statistical views with ability to export this data to MS Excel ®. In all formats the data has been validated to be correct.

DATE

1 January, 2009

Approved signatory typed name and position:

Dov Bruker



Fourier Systems (1989) Ltd.

CEO

Fourier Systems

Rosh Ha'ayin, Israel

Tel: +972-3-901-4849

[www.fouriersystems.com](http://www.fouriersystems.com)

Manufactured under ISO 9001:2000

